

Filtek™ Z250 Universal Restorative

Reduced shrinkage, proven performance and excellent esthetics.

micro-hybrid



The lowest polymerization shrinkage (see chart) is only the beginning. With its high performance, strength and easy handling, this micro-hybrid restorative fills a variety of needs—with beautiful results.

Filtek™ Z250 Universal Restorative offers:

Lowest Polymerization Shrinkage (see chart) – Less stress on the tooth, less sensitivity potential for the patient.

Speed and Convenience – Easy-to-use, direct-delivery capsule. Cures in just 20 seconds.* Easy to finish to a smooth surface.

Excellent Handling – Exceptional resistance to slumping.


Outstanding Versatility – Useful for any area of the mouth.

Excellent Esthetics – Available in 15 shades to meet your needs.

Filtek Z250 universal restorative gives you strong, long-lasting, natural-looking results.

KEY FEATURES

- LOW SHRINKAGE
- EXCELLENT ESTHETICS
- ALL-AROUND VERSATILITY



Lowest shrinkage of major competitors. 3M ESPE R&D conducted a shrinkage test on a variety of market-leading restoratives. Of those tested Filtek™ Z250 Universal Restorative demonstrated the least amount of polymerization shrinkage.



Filtek™ Z250 Universal Restorative

Ordering Information

3410 3M™ ESPE™ Single Bond Dental Adhesive System



No. 3410 Includes:

- 2 - 3ml phosphoric etchant syringes
- 6ml adhesive
- Brush handle, brush tips, etchant syringe tips, disposable mixing wells, technique guides

- Filtek Z250 Universal Restorative and the popular Single Bond dental adhesive work together as a system to give you excellent performing restorations.
- Independent studies give Single Bond dental adhesive excellent ratings for ease of use and high bond strength properties.

Single Bond Dental Adhesive Refills

3411	Single Bond Adhesive refill, 6ml
7523	3M™ ESPE™ Scotchbond™ Phosphoric Etchant, 2-3ml syringes

6020 Filtek Z250 Universal Restorative System Syringe Introductory Kit



No. 6020 Introductory Kit Includes:

- 1 - 4g syringe each of Filtek Z250 universal restorative in shades A2, A3, B2 and C3
- 1 - 3ml vial of Single Bond dental adhesive
- 1 - 3ml syringe of Scotchbond etchant and syringe tips
- Accessories

6021 Filtek Z250 Universal Restorative System Capsule Introductory Kit



No. 6021 Introductory Kit Includes:

- 10 - .20g capsules each of Filtek Z250 universal restorative in shades A2, A3, B0.5, B2 and C3
- 1 - 3ml vial of Single Bond dental adhesive
- 1 - 3ml syringe of Scotchbond etchant and syringe tips
- Accessories

Filtek Z250 Universal Restorative Refills

Shades	Item Number	
	Syringe	Capsule
A1	6020A1	6021A1
A2	6020A2	6021A2
A3	6020A3	6021A3
A3.5	6020A3.5	6021A3.5
A4	6020A4	6021A4
B0.5	6020B0.5	6021B0.5
B1	6020B1	6021B1
B2	6020B2	6021B2
B3	6020B3	6021B3
C2	6020C2	6021C2
C3	6020C3	6021C3
C4	6020C4	6021C4
D3	6020D3	6021D3
Incisal	6020I	6021I
Universal Dentin	6020UD	6021UD

Indications for use

	Filtek Z250 Universal Restorative
Class I	•
Class II	•
Class III	•
Class IV	•
Class V	•
Veneer	•
Pedo	•

3M, ESPE, Filtek and Scotchbond are trademarks of 3M or 3M ESPE AG. Herculite XRV and Prodigy are trademarks of Kerr Corporation. Charisma is a registered trademark of Heraeus Kulzer, Inc. Tetric Ceram is a trademark of Ivoclar Vivadent. TPH is a trademark of Dentsply.

*Shades UD, C4, and B0.5 require a 30-second cure for a 2mm increment.

3M ESPE Technical Hotline: 1-800-634-2249

3M ESPE Web site: www.3MESPE.com

3M ESPE

Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
USA

Dental Products

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1 800 265-1840 ext. 6229



Minimum 10%
Post-Consumer Fiber

Printed in U.S.A.
© 3M 1999-2002
70-2009-2191-7

Ordering